

Product datasheet for **HP206459**

PTP alpha (PTPRA) Human qPCR Primer Pair (NM_002836)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5786
Forward Sequence:	GGTGTCTGTAGAGGATGTGACTG
Reverse Sequence:	AGTGATGAGGCGCTGTGGCTTT
ACCN:	<u>NM_002836</u> , <u>NM_002836.1</u> , <u>NM_002836.2</u> , <u>NM_002836.3</u> , <u>BC027308</u> , <u>BC027308.1</u> , <u>BI597187</u> , <u>BQ213974</u> , <u>BX571753</u>
UniProt ID:	<u>P18433</u>
Synonyms:	HEPTP; HLPR; HPTPA; HPTPalpha; LRP; PTPA; PTPRL2; R-PTP-alpha; RPTPA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)