

Product datasheet for **HP206475**

PCYT2 Human qPCR Primer Pair (NM_002861)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5833
Forward Sequence:	GGGCTTACTACTTTGACCAGGAG
Reverse Sequence:	GGCTCCAATCACCACCTTCTGAC
ACCN:	NM_002861 , NM_002861.1 , NM_002861.3 , NM_002861.4 , BC010075 , BC010075.2 , BC000351 , BC073815 , BX419403 , BX459568 , NM_002861.5
UniProt ID:	Q99447
Synonyms:	ET; SPG82
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)