

## Product datasheet for **HP206497**

### **RANGAP1 Human qPCR Primer Pair (NM\_002883)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	5905
<b>Forward Sequence:</b>	GGTGGTCTCTGCCTTCTAAAG
<b>Reverse Sequence:</b>	CGAGGAGTTGAAAGCCTTCTGC
<b>ACCN:</b>	<a href="#">NM_002883</a> , <a href="#">NM_002883.1</a> , <a href="#">NM_002883.2</a> , <a href="#">NM_002883.3</a> , <a href="#">BC014044</a> , <a href="#">BC014044.2</a> , <a href="#">BC004990</a> , <a href="#">BC019676</a> , <a href="#">BC041396</a> , <a href="#">BC048990</a> , <a href="#">BM551472</a> , <a href="#">NM_002883.4</a>
<b>UniProt ID:</b>	<a href="#">P46060</a>
<b>Synonyms:</b>	Fug1; RANGAP; SD
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)