

Product datasheet for **HP206498**

RAP1GAP Human qPCR Primer Pair (NM_002885)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5909
Forward Sequence:	CCGACATGATCGCGTCCAACCTT
Reverse Sequence:	GGGTCCAAGAAGGGCACATCA
ACCN:	<u>NM_002885</u> , <u>NM_002885.1</u> , <u>NM_002885.2</u> , <u>BC054490</u> , <u>BC035030</u> , <u>BM739021</u> , <u>BP200029</u> , <u>BP283134</u> , <u>NM_002885.4</u>
UniProt ID:	<u>P47736</u>
Synonyms:	RAP1GA1; RAP1GAP1; RAP1GAPII; RAPGAP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)