

Product datasheet for **HP206499**

RAP2B Human qPCR Primer Pair (NM_002886)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5912
Forward Sequence:	CGCAAGGAGATTGAGGTGGACT
Reverse Sequence:	TTGACGAGGCTGTAGACCAGGA
ACCN:	<u>NM_002886</u> , <u>NM_002886.1</u> , <u>NM_002886.3</u> , <u>BC012362</u> , <u>BC012362.1</u> , <u>BF965860</u> , <u>BI551390</u> , <u>BQ378964</u> , <u>BU507815</u> , <u>BU860799</u> , <u>BX116683</u> , <u>BX355341</u> , <u>NM_002886.4</u>
UniProt ID:	<u>P61225</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)