

Product datasheet for **HP206503**

RASA1 Human qPCR Primer Pair (NM_002890)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5921
Forward Sequence:	GGGACATCCAATAAACGCCTTCG
Reverse Sequence:	TTTGCTACTTGGACACTATTCAGG
ACCN:	<u>NM_002890</u> , <u>NM_002890.1</u> , <u>NM_002890.2</u> , <u>BC033015</u> , <u>BC033015.2</u> , <u>BC020761</u> , <u>BC054891</u> , <u>BC070063</u> , <u>BC107870</u> , <u>BC131816</u> , <u>NM_002890.3</u>
UniProt ID:	<u>P20936</u>
Synonyms:	CM-AVM; CMAVM; CMAVM1; GAP; p120; p120GAP; p120RASGAP; PKWS; RASA; RASGAP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)