

## Product datasheet for **HP206507**

### CtIP (RBBP8) Human qPCR Primer Pair (NM\_002894)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5932
Forward Sequence:	TGGCAGACAGTTTCTCCCAAGC
Reverse Sequence:	GGCTCCACAAACGCTTTCTGCT
ACCN:	<a href="#">NM_002894</a> , <a href="#">NM_002894.2</a> , <a href="#">BC001170</a> , <a href="#">BC030590</a> , <a href="#">BG723261</a> , <a href="#">BU187672</a> , <a href="#">BX648221</a> , <a href="#">NM_002894.3</a>
UniProt ID:	<a href="#">Q99708</a>
Synonyms:	COM1; CTIP; JWDS; RIM; SAE2; SCKL2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)