

Product datasheet for **HP206524**

UPF1 Human qPCR Primer Pair (NM_002911)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5976
Forward Sequence:	AACGAGCACCAAGGCATTGGCT
Reverse Sequence:	GGCTGCTTTGATAGTGCCTTCG
ACCN:	<u>NM_002911</u> , <u>NM_002911.1</u> , <u>NM_002911.2</u> , <u>NM_002911.3</u> , <u>BC039817</u> , <u>BC039817.1</u> , <u>BC034291</u> , <u>BE243084</u> , <u>BM969672</u> , <u>BQ950083</u> , <u>BX331580</u> , <u>BX364728</u> , <u>NM_002911.4</u>
UniProt ID:	<u>Q92900</u>
Synonyms:	HUPF1; NORF1; pNORF1; RENT1; smg-2; UTF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)