

Product datasheet for **HP206539**

RGS16 Human qPCR Primer Pair (NM_002928)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6004
Forward Sequence:	GAAGATCCGATCAGCTACCAAGC
Reverse Sequence:	TGCAGGTTTCATCCTCGTCAGCT
ACCN:	NM_002928 , NM_002928.1 , NM_002928.2 , NM_002928.3 , BC006243 , BC006243.2 , BT006638
UniProt ID:	O15492
Synonyms:	A28-RGS14; A28-RGS14P; RGS-R
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)