

Product datasheet for **HP206550**

ROBO1 Human qPCR Primer Pair (NM_002941)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6091
Forward Sequence:	AGTGAGCCTCAGTTCATCCAGC
Reverse Sequence:	GCTCCAATACCTGCTATGAAGGC
ACCN:	NM_002941 , NM_002941.1 , NM_002941.2 , NM_002941.3 , BC157861 , BC115022 , BC001969 , BC029676 , BC112336 , BC115020 , BC171855 , BP228659 , BX538319 , NM_002941.4
UniProt ID:	Q9Y6N7
Synonyms:	DUTT1; SAX3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)