

Product datasheet for **HP206567**

S100A4 Human qPCR Primer Pair (NM_002961)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6275
Forward Sequence:	CAGAACTAAAGGAGCTGCTGACC
Reverse Sequence:	CTTGGAAGTCCACCTCGTTGTC
ACCN:	<u>NM_002961</u> , <u>NM_002961.1</u> , <u>NM_002961.2</u> , <u>BC000838</u> , <u>BC000838.1</u> , <u>BC016300</u> , <u>BC016300.1</u> , <u>BC024974</u> , <u>NM_002961.3</u>
UniProt ID:	<u>P26447</u>
Synonyms:	18A2; 42A; CAPL; FSP1; MTS1; P9KA; PEL98
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)