

Product datasheet for **HP206569**

Psoriasis (S100A7) Human qPCR Primer Pair (NM_002963)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6278
Forward Sequence:	AGAAGCCAAGCCTGCTGACGAT
Reverse Sequence:	GTCCTTTTCTCAAAGACATCGGC
ACCN:	NM_002963 , NM_002963.1 , NM_002963.2 , NM_002963.3 , BC034687 , BC034687.1 , BF090388 , BX476368 , NM_002963.4
UniProt ID:	P31151
Synonyms:	PSOR1; S100A7c
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)