

Product datasheet for **HP206571**

Calprotectin (S100A9) Human qPCR Primer Pair (NM_002965)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6280
Forward Sequence:	GCACCCAGACACCCTGAACCA
Reverse Sequence:	TGTGTCCAGGTCCTCCATGATG
ACCN:	NM_002965 , NM_002965.1 , NM_002965.2 , NM_002965.3 , BC047681 , BC047681.2 , NM_002965.4
UniProt ID:	P06702
Synonyms:	60B8AG; CAGB; CFAG; CGLB; L1AG; LIAG; MAC387; MIF; MRP14; NIF; P14
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)