

## Product datasheet for **HP206572**

### S100A10 Human qPCR Primer Pair (NM\_002966)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6281
Forward Sequence:	AACAAAGGAGGACCTGAGAGTAC
Reverse Sequence:	CTTGCCATCTCTACACTGGTCC
ACCN:	<u><a href="#">NM_002966</a></u> , <u><a href="#">NM_002966.1</a></u> , <u><a href="#">NM_002966.2</a></u> , <u><a href="#">BC015973</a></u> , <u><a href="#">BC015973.1</a></u> , <u><a href="#">BC105786</a></u> , <u><a href="#">BE378444</a></u> , <u><a href="#">BI827113</a></u> , <u><a href="#">BQ050753</a></u> , <u><a href="#">NM_002966.3</a></u>
UniProt ID:	<u><a href="#">P60903</a></u>
Synonyms:	42C; ANX2L; ANX2LG; CAL1L; Ca[1]; CLP11; GP11; p10; P11
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)