

Product datasheet for **HP206573**

SAFB Human qPCR Primer Pair (NM_002967)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6294
Forward Sequence:	TAACCACCTGCACAAGACGGAG
Reverse Sequence:	GTCACTGTTGCTCGACTTCTCC
ACCN:	NM_002967 , NM_002967.1 , NM_002967.2 , NM_002967.3 , BC020309 , BC032055 , BC126219 , BC143937 , BC143938 , BC143939
UniProt ID:	Q15424
Synonyms:	HAP; HET; SAB-B1; SAF-B; SAF-B1; SAFB1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)