

Product datasheet for **HP206574**

SALL1 Human qPCR Primer Pair (NM_002968)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6299
Forward Sequence:	GCTTGCACTATTTGTGGAAGAGC
Reverse Sequence:	GAACCTGACGGGATTGCCTCCT
ACCN:	<u>NM_002968</u> , <u>NM_002968.1</u> , <u>NM_002968.2</u> , <u>BC113881</u> , <u>BC113906</u> , <u>BC131548</u> , <u>BC167772</u> , <u>BP229682</u> , <u>NM_002968.3</u>
UniProt ID:	<u>Q9NSC2</u>
Synonyms:	HEL-S-89; HSAL1; Sal-1; TBS; ZNF794
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)