

Product datasheet for **HP206575**

SAPK3 (MAPK12) Human qPCR Primer Pair (NM_002969)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6300
Forward Sequence:	CAGTTCCTCGTGACCAGATGC
Reverse Sequence:	CACAGTCTTCGTTACAGCCAG
ACCN:	NM_002969 , NM_002969.1 , NM_002969.2 , NM_002969.3 , NM_002969.4 , NM_002969.5 , BC015741 , BC015741.1 , BX640648 , NM_002969.6
UniProt ID:	P53778
Synonyms:	ERK-6; ERK3; ERK6; MAPK 12; P38GAMMA; PRKM12; SAPK-3; SAPK3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)