

## Product datasheet for **HP206587**

### MCP1 (CCL2) Human qPCR Primer Pair (NM\_002982)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6347
Forward Sequence:	AGAATCACCAGCAGCAAGTGTC
Reverse Sequence:	TCCTGAACCCACTTCTGCTTG
ACCN:	<a href="#">NM_002982</a> , <a href="#">NM_002982.1</a> , <a href="#">NM_002982.2</a> , <a href="#">NM_002982.3</a> , <a href="#">BC009716</a> , <a href="#">BC009716.1</a> , <a href="#">BG530064</a> , <a href="#">BT007329</a> , <a href="#">NM_002982.4</a>
UniProt ID:	<a href="#">P13500</a>
Synonyms:	GDCF-2; HC11; HSMCR30; MCAF; MCP-1; MCP1; SCYA2; SMC-CF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)