

Product datasheet for **HP206590**

RANTES (CCL5) Human qPCR Primer Pair (NM_002985)

Product data:

| | |
|--------------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 6352 |
| Forward Sequence: | CCTGCTGCTTTGCCTACATTGC |
| Reverse Sequence: | ACACACTTGCGGTTCTTTCCGG |
| ACCN: | BC003528 , BC007399 , BC009519 , NM_002985 , NM_002985.1 , NM_002985.2 , BC008600 , BG272739 , NM_002985.3 |
| UniProt ID: | P13501 |
| Synonyms: | D17S136E; eoCP; RANTES; SCYA5; SIS-delta; SISd; TCP228 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)