

## Product datasheet for **HP206610**

### SELENOW Human qPCR Primer Pair (NM\_003009)

#### Product data:

|                   |                                                                                                                                                                                                                                                                                                                                               |
|-------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:     | qPCR Primer Pairs                                                                                                                                                                                                                                                                                                                             |
| Gene ID:          | 6415                                                                                                                                                                                                                                                                                                                                          |
| Forward Sequence: | GGTTCTTTGAAGTGATGGTAGCC                                                                                                                                                                                                                                                                                                                       |
| Reverse Sequence: | TGAGCCAAGGCGGCTTTGATGG                                                                                                                                                                                                                                                                                                                        |
| ACCN:             | <a href="#">NM_003009</a> , <a href="#">NM_003009.1</a> , <a href="#">NM_003009.2</a> , <a href="#">BC039597</a> , <a href="#">BC039597.1</a> , <a href="#">BC000581</a> , <a href="#">BC032546</a> , <a href="#">BC047893</a> , <a href="#">BF032577</a> , <a href="#">BF973456</a> , <a href="#">BU185196</a> , <a href="#">NM_003009.4</a> |
| UniProt ID:       | <a href="#">P63302</a>                                                                                                                                                                                                                                                                                                                        |
| Synonyms:         | selW; SEPW1                                                                                                                                                                                                                                                                                                                                   |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.                                                                                                                                   |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.                                                                                                                                                                                                                                                                                                                               |
| Stability:        | The primer mix is stable for one year from date of shipping.                                                                                                                                                                                                                                                                                  |



[View online »](#)