

Product datasheet for **HP206613**

SFRP2 Human qPCR Primer Pair (NM_003013)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6423
Forward Sequence:	CTCCAAAGGTATGTGAAGCCTGC
Reverse Sequence:	CCAGGATGATTTTGGTATCTCGG
ACCN:	NM_003013 , NM_003013.1 , NM_003013.2 , BC008666 , BC008666.1 , BM906663 , BM924815 , BQ130455 , NM_003013.3
UniProt ID:	Q96HF1
Synonyms:	FRP-2; SARP1; SDF-5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)