

Product datasheet for **HP206614**

SFRP4 Human qPCR Primer Pair (NM_003014)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6424
Forward Sequence:	CTATGACCGTGGCGTGTGCATT
Reverse Sequence:	GCTTAGGGGTTTACAGTCAACATC
ACCN:	<u>NM_003014</u> , <u>NM_003014.2</u> , <u>NM_003014.3</u> , <u>BC032828</u> , <u>BC047684</u> , <u>BC058911</u> , <u>BT019679</u>
UniProt ID:	<u>Q6FHJ7</u>
Synonyms:	FRP-4; FRPHE; FRZB-2; PYL; sFRP-4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)