

Product datasheet for **HP206616**

SC35 (SRSF2) Human qPCR Primer Pair (NM_003016)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6427
Forward Sequence:	GCGGTGGCTACGGACGCCG
Reverse Sequence:	TTCGAGCGGCTGTAGCGAGATC
ACCN:	<u>NM_003016</u> , <u>NM_003016.1</u> , <u>NM_003016.2</u> , <u>NM_003016.3</u> , <u>NM_003016.4</u> , <u>BC070086</u> , <u>BC000339</u> , <u>BC001303</u> , <u>BC063570</u> , <u>BC066958</u> , <u>BT007250</u>
UniProt ID:	<u>Q01130</u>
Synonyms:	PR264; SC-35; SC35; SFRS2; SFRS2A; SRp30b
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)