

Product datasheet for **HP206618**

SFTPC Human qPCR Primer Pair (NM_003018)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6440
Forward Sequence:	GTCCTCATCGTCGGTGATTG
Reverse Sequence:	AGAAGGTGGCAGTGGTAACCAG
ACCN:	<u>NM_003018</u> , <u>NM_003018.1</u> , <u>NM_003018.2</u> , <u>NM_003018.3</u> , <u>BC005913</u> , <u>BC005913.1</u> , <u>BC105594</u> , <u>BF433929</u> , <u>BG535994</u> , <u>BG569816</u> , <u>BI490051</u> , <u>BI524154</u> , <u>BM993053</u> , <u>BQ027183</u>
UniProt ID:	<u>P11686</u>
Synonyms:	BRICD6; PSP-C; SFTP2; SMDP2; SP-C; SP5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)