

Product datasheet for **HP206620**

Secretogranin V (SCG5) Human qPCR Primer Pair (NM_003020)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6447
Forward Sequence:	CGGAACATGACTATCCAGGCTTG
Reverse Sequence:	CATTATCCAGTCTCTGTCCTTGTA
ACCN:	NM_003020 , NM_003020.1 , NM_003020.2 , NM_003020.3 , NM_003020.4 , BC005349 , BC005349.1 , BC093053 , BP308267 , BP377248
UniProt ID:	P05408
Synonyms:	7B2; P7B2; SGNE1; SgV
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)