

Product datasheet for **HP206624**

SH3GL2 Human qPCR Primer Pair (NM_003026)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6456
Forward Sequence:	ACCACAAGCAGGCAGTCCAGAT
Reverse Sequence:	GGAAACTCCAGGCTCATTCTGTG
ACCN:	<u>NM_003026</u> , <u>NM_003026.1</u> , <u>NM_003026.2</u> , <u>NM_003026.3</u> , <u>BC032825</u> , <u>BT019682</u> , <u>BX538074</u> , <u>NM_003026.5</u>
UniProt ID:	<u>Q99962</u>
Synonyms:	CNSA2; EEN-B1; SH3D2A; SH3P4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)