

Product datasheet for **HP206643**

SLC10A1 Human qPCR Primer Pair (NM_003049)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6554
Forward Sequence:	GCTCTCTTCTGCCTCAATGGAC
Reverse Sequence:	AGTGGTCCAATGACTTCAGGTGG
ACCN:	<u>NM_003049</u> , <u>NM_003049.1</u> , <u>NM_003049.2</u> , <u>NM_003049.3</u> , <u>BC069799</u> , <u>BC069822</u> , <u>BC074724</u> , <u>BC126298</u> , <u>BC136355</u> , <u>NM_003049.4</u>
UniProt ID:	<u>Q14973</u>
Synonyms:	FHCA2; NTCP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)