

## Product datasheet for **HP206655**

### SLIT3 Human qPCR Primer Pair (NM\_003062)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6586
Forward Sequence:	CCAGATGACTGTGAGGACAACG
Reverse Sequence:	CACAGTGGTCAATCACCTCGTC
ACCN:	<a href="#">NM_003062</a> , <a href="#">NM_003062.1</a> , <a href="#">NM_003062.2</a> , <a href="#">NM_003062.3</a> , <a href="#">BC062365</a> , <a href="#">BC098388</a> , <a href="#">BC142665</a> , <a href="#">BC142717</a> , <a href="#">BC146759</a> , <a href="#">NM_003062.4</a>
UniProt ID:	<a href="#">O75094</a>
Synonyms:	MEGF5; SLIL2; Slit-3; SLIT1; slit2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)