

Product datasheet for **HP206660**

SMARCA2 Human qPCR Primer Pair (NM_003070)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6595
Forward Sequence:	GAAGAGTCCAGTAGGCAGGAAAC
Reverse Sequence:	CATCCACGTCTTGCTTAGCTGTC
ACCN:	<u>NM_003070</u> , <u>NM_003070.1</u> , <u>NM_003070.2</u> , <u>NM_003070.3</u> , <u>NM_003070.4</u> , <u>BC040029</u> , <u>BC068252</u> , <u>BC156185</u> , <u>BI758293</u> , <u>BM921013</u> , <u>NM_003070.5</u>
UniProt ID:	<u>P51531</u>
Synonyms:	BAF190; BIS; BRM; hBRM; hSNF2a; NCBRS; SNF2; SNF2L2; SNF2LA; Sth1p; SWI2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)