

Product datasheet for **HP206689**

SorLA (SORL1) Human qPCR Primer Pair (NM_003105)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6653
Forward Sequence:	GAACACCTGTCTTCGCAACCAG
Reverse Sequence:	TGTCCAGGTCACAGATGGTGGT
ACCN:	<u>NM_003105</u> , <u>NM_003105.1</u> , <u>NM_003105.2</u> , <u>NM_003105.3</u> , <u>NM_003105.4</u> , <u>NM_003105.5</u> , <u>BC032903</u> , <u>BC040643</u> , <u>BC137171</u> , <u>BQ888333</u> , <u>BX365721</u> , <u>NM_003105.6</u>
UniProt ID:	<u>Q92673</u>
Synonyms:	C11orf32; gp250; LR11; LRP9; SORLA; SorLA-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)