

Product datasheet for **HP206690**

SOX2 Human qPCR Primer Pair (NM_003106)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6657
Forward Sequence:	GCTACAGCATGATGCAGGACCA
Reverse Sequence:	TCTGCGAGCTGGTCATGGAGTT
ACCN:	<u>BC013923</u> , <u>NM_003106</u> , <u>NM_003106.1</u> , <u>NM_003106.2</u> , <u>NM_003106.3</u> , <u>BC013923.2</u> , <u>BF305585</u> , <u>BM668019</u> , <u>BM722297</u> , <u>NM_003106.4</u>
UniProt ID:	<u>P48431</u>
Synonyms:	ANOP3; MCOPS3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)