

Product datasheet for **HP206695**

SPAG1 Human qPCR Primer Pair (NM_003114)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6674
Forward Sequence:	GCAGATAGACTGTGGACTCCAG
Reverse Sequence:	GCAGGAATAGGTGACAGCTTCTC
ACCN:	<u>NM_003114</u> , <u>NM_003114.1</u> , <u>NM_003114.2</u> , <u>NM_003114.3</u> , <u>BC055091</u> , <u>BM557203</u> , <u>BM981630</u> , <u>BX480304</u> , <u>NM_003114.5</u>
UniProt ID:	<u>Q07617</u>
Synonyms:	CILD28; CT140; DNAAF13; HEL-S-268; HSD-3.8; SP75; TPIS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)