

## Product datasheet for **HP206700**

### SPG7 Human qPCR Primer Pair (NM\_003119)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6687
Forward Sequence:	AGCAGACGCTCAACCAGCTTCT
Reverse Sequence:	CCTCATCAGAGCACCGTCCAAA
ACCN:	<a href="#">NM_003119</a> , <a href="#">NM_003119.1</a> , <a href="#">NM_003119.2</a> , <a href="#">NM_003119.3</a> , <a href="#">BC036104</a> , <a href="#">BC036104.1</a> , <a href="#">BC007692</a> , <a href="#">BC015411</a> , <a href="#">BC035929</a> , <a href="#">BC052587</a> , <a href="#">BI827568</a> , <a href="#">NM_003119.4</a>
UniProt ID:	<a href="#">Q9UQ90</a>
Synonyms:	CAR; CMAR; PGN; SPG5C
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)