

Product datasheet for **HP206722**

SSBP1 Human qPCR Primer Pair (NM_003143)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6742
Forward Sequence:	CTGTCTTGAGACAGGTGGAAGG
Reverse Sequence:	CTGTGCCATGTTGTCTTTTACTG
ACCN:	NM_003143 , NM_003143.2 , BC000895 , BC008402 , BC093054 , BF131408 , BG030060 , BP308841 , BQ959702 , BU532904 , NM_003143.3
UniProt ID:	Q04837
Synonyms:	Mt-SSB; mtSSB; SOSS-B1; SSBP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)