

Product datasheet for **HP206754**

ADAM17 Human qPCR Primer Pair (NM_003183)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6868
Forward Sequence:	AACAGCGACTGCACGTTGAAGG
Reverse Sequence:	CTGTGCAGTAGGACACGCCTTT
ACCN:	NM_003183 , NM_003183.4 , BC062687 , BC136783 , BC146658 , BM725368 , BQ186514 , BU734977 , BX339672
UniProt ID:	P78536
Synonyms:	ADAM18; CD156B; CSVP; NISBD; NISBD1; TACE
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)