

Product datasheet for **HP206755**

TAF2 Human qPCR Primer Pair (NM_003184)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6873
Forward Sequence:	CATCTCGGCTTGCACTCACTGA
Reverse Sequence:	AGAGTGACTTCATGGCTGGTGG
ACCN:	NM_003184 , NM_003184.1 , NM_003184.2 , NM_003184.3 , BC008473 , BC018010 , BC035673 , BC037262 , BC047732 , BC064830 , BC146379
UniProt ID:	Q6P1X5
Synonyms:	CIF150; MRT40; TAF2B; TAFII150
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)