

Product datasheet for **HP206762**

TATA binding protein (TBP) Human qPCR Primer Pair (NM_003194)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6908
Forward Sequence:	TGTATCCACAGTGAATCTTGTTG
Reverse Sequence:	GGTTCGTGGCTCTTATCCTC
ACCN:	<u>NM_003194</u> , <u>NM_003194.1</u> , <u>NM_003194.2</u> , <u>NM_003194.3</u> , <u>BC110341</u> , <u>BC110341.1</u> , <u>BC109053</u> , <u>BC109054</u> , <u>BM720195</u> , <u>BT019657</u> , <u>NM_003194.5</u>
UniProt ID:	<u>P20226</u>
Synonyms:	GTF2D; GTF2D1; HDL4; SCA17; TFIID
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)