

Product datasheet for **HP206763**

TFIIS (TCEA2) Human qPCR Primer Pair (NM_003195)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6919
Forward Sequence:	AACAAGTGCCGCGAGATGCTGA
Reverse Sequence:	GATACGACTCCGTACACGGTTC
ACCN:	<u>NM_003195</u> , <u>NM_003195.1</u> , <u>NM_003195.2</u> , <u>NM_003195.3</u> , <u>NM_003195.4</u> , <u>NM_003195.5</u> , <u>BC018896</u> , <u>BC018896.2</u> , <u>BC014211</u> , <u>BC031877</u> , <u>BC050623</u> , <u>BC050624</u> , <u>BC056407</u> , <u>BC098585</u> , <u>BI668232</u> , <u>BI756937</u> , <u>NM_003195.6</u>
UniProt ID:	<u>Q15560</u>
Synonyms:	TFIIS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)