

Product datasheet for **HP206769**

TCF7 Human qPCR Primer Pair (NM_003202)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6932
Forward Sequence:	CTGACCTCTCTGGCTTCTACTC
Reverse Sequence:	CAGAACCTAGCATCAAGGATGGG
ACCN:	NM_003202 , NM_003202.1 , NM_003202.2 , NM_003202.3 , NM_003202.4 , BC048769 , BC048769.1 , BC072023 , NM_003202.5
UniProt ID:	P36402
Synonyms:	TCF-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)