

Product datasheet for **HP206775**

TEAD4 Human qPCR Primer Pair (NM_003213)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7004
Forward Sequence:	GAAGGTCTGCTCTTTCGGCAAG
Reverse Sequence:	GAGGTGCTTGAGCTTGTGGATG
ACCN:	<u>NM_003213</u> , <u>NM_003213.1</u> , <u>NM_003213.2</u> , <u>NM_003213.3</u> , <u>BC015497</u> , <u>BC015497.1</u> , <u>BM464237</u> , <u>BM555991</u> , <u>BX640921</u> , <u>NM_003213.4</u>
UniProt ID:	<u>Q15561</u>
Synonyms:	EFTR-2; hRTEF-1B; RTEF1; TCF13L1; TEF-3; TEF3; TEFR-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)