

Product datasheet for **HP206776**

TEAD3 Human qPCR Primer Pair (NM_003214)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 7005 |
| Forward Sequence: | AGGCAGTAGATGTGCGCCAGAT |
| Reverse Sequence: | TCCTGGATGGTGTGTTGAGGT |
| ACCN: | <u>NM_003214</u> , <u>NM_003214.1</u> , <u>NM_003214.2</u> , <u>NM_003214.3</u> , <u>BC091488</u> , <u>BC091488.1</u> , <u>BC010053</u> , <u>BC027877</u> , <u>BQ719548</u> , <u>NM_003214.4</u> |
| UniProt ID: | <u>Q99594</u> |
| Synonyms: | DTEF-1; ETEF-1; TEAD-3; TEAD5; TEF-5; TEF5 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)