

Product datasheet for **HP206784**

ARFRP1 Human qPCR Primer Pair (NM_003224)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10139
Forward Sequence:	GTGCCTCTCAATCCCTGACATC
Reverse Sequence:	ATTCCGCACGACACACTTCACC
ACCN:	<u>NM_003224</u> , <u>NM_003224.1</u> , <u>NM_003224.2</u> , <u>NM_003224.3</u> , <u>NM_003224.4</u> , <u>NM_003224.5</u> , <u>BC092480</u> , <u>BC092480.1</u> , <u>BC021023</u> , <u>BC093967</u> , <u>BC143447</u> , <u>BC143448</u> , <u>BC143467</u> , <u>BC143472</u> , <u>BC143477</u> , <u>BI524032</u> , <u>BI711667</u> , <u>BI819763</u> , <u>BM546222</u> , <u>BM549361</u> , <u>BP384879</u> , <u>BT019624</u> , <u>BT019663</u> , <u>BU526941</u> , <u>BX342690</u> , <u>BX342805</u> , <u>BX463457</u> , <u>NM_003224.6</u>
UniProt ID:	<u>Q13795</u>
Synonyms:	ARL18; ARP; Arp1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)