

Product datasheet for **HP206785**

TFF1 Human qPCR Primer Pair (NM_003225)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 7031 |
| Forward Sequence: | CCAGTGTGCAAATAAGGGCTGC |
| Reverse Sequence: | AGGCAGATCCCTGCAGAAGTGT |
| ACCN: | <u>NM_003225</u> , <u>NM_003225.1</u> , <u>NM_003225.2</u> , <u>BC032811</u> , <u>BC032811.1</u> , <u>BQ083523</u> , <u>BU539534</u> , <u>NM_003225.3</u> |
| UniProt ID: | <u>P04155</u> |
| Synonyms: | BCEI; D21S21; HP1.A; HPS2; pNR-2; pS2 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)