

## Product datasheet for **HP206790**

### TGF beta 2 (TGFB2) Human qPCR Primer Pair (NM\_003238)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7042
Forward Sequence:	AAGAAGCGTGCTTTGGATGCGG
Reverse Sequence:	ATGCTCCAGCACAGAAGTTGGC
ACCN:	<a href="#">NM_003238</a> , <a href="#">NM_003238.1</a> , <a href="#">NM_003238.2</a> , <a href="#">NM_003238.3</a> , <a href="#">NM_003238.4</a> , <a href="#">BC096235</a> , <a href="#">BC096236</a> , <a href="#">BC099635</a> , <a href="#">BC099636</a> , <a href="#">BX098793</a> , <a href="#">NM_003238.6</a>
UniProt ID:	<a href="#">P61812</a>
Synonyms:	G-TSF; LDS4; TGF-beta2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)