

Product datasheet for **HP206804**

TIMP1 Human qPCR Primer Pair (NM_003254)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 7076 |
| Forward Sequence: | GGAGAGTGTCTGCGGATACTTC |
| Reverse Sequence: | GCAGGTAGTGATGTGCAAGAGTC |
| ACCN: | <u>NM_003254</u> , <u>NM_003254.1</u> , <u>NM_003254.2</u> , <u>BC000866</u> , <u>BC000866.1</u> , <u>BC007097</u> , <u>BQ181804</u> , <u>BU857950</u> |
| UniProt ID: | <u>P01033</u> |
| Synonyms: | CLGI; EPA; EPO; HCI; TIMP; TIMP-1 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)