

## Product datasheet for **HP206811**

### SEC62 Human qPCR Primer Pair (NM\_003262)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7095
Forward Sequence:	TCAGTGTGGGTGCAGGCTGTTT
Reverse Sequence:	ACCAAAAGTGGTGCCTTCCTCC
ACCN:	<a href="#">NM_003262</a> , <a href="#">NM_003262.1</a> , <a href="#">NM_003262.2</a> , <a href="#">NM_003262.3</a> , <a href="#">BC012035</a> , <a href="#">BC012035.1</a> , <a href="#">BC003677</a> , <a href="#">BC007040</a> , <a href="#">BC089399</a> , <a href="#">BC105971</a> , <a href="#">BX475495</a> , <a href="#">BX537668</a> , <a href="#">BX648539</a> , <a href="#">NM_003262.4</a>
UniProt ID:	<a href="#">Q99442</a>
Synonyms:	Dtrp1; HTP1; TLOC1; TP-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)