

Product datasheet for **HP206818**

TSPAN6 Human qPCR Primer Pair (NM_003270)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7105
Forward Sequence:	CGTGCTCATTGCTACTGGTACC
Reverse Sequence:	CAAATCCTACGATGGCAGCGAC
ACCN:	<u>NM_003270</u> , <u>NM_003270.1</u> , <u>NM_003270.2</u> , <u>NM_003270.3</u> , <u>BC012389</u> , <u>BC012389.1</u> , <u>BI458635</u> , <u>BT006977</u> , <u>NM_003270.4</u>
UniProt ID:	<u>O43657</u>
Synonyms:	T245; TM4SF6; TSPAN-6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)