

Product datasheet for **HP206828**

Troponin T1 (TNNT1) Human qPCR Primer Pair (NM_003283)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7138
Forward Sequence:	AACGCGAACGTCAGGCTAAGCT
Reverse Sequence:	CTTGACCAGGTAGCCGCCAAAA
ACCN:	NM_003283 , NM_003283.1 , NM_003283.2 , NM_003283.3 , NM_003283.4 , NM_003283.5 , BC107798 , BC010963 , BC007074 , BC022086 , BC034143 , BT019630 , NM_003283.6
UniProt ID:	P13805
Synonyms:	ANM; NEM5; STNT; TNT; TNTS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)