

Product datasheet for **HP206839**

TCTP (TPT1) Human qPCR Primer Pair (NM_003295)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7178
Forward Sequence:	TTGGACTACCGTGAGGATGGTG
Reverse Sequence:	CAGAAGCCAGTTATGATGACAGG
ACCN:	BC022436 , NM_003295 , NM_003295.1 , NM_003295.2 , NM_003295.3 , BC003352 , BC003352.2 , BC007842 , BC012431 , BC040008 , BC052333 , BP215975 , BC022436.1 , NM_003295.4
UniProt ID:	P13693
Synonyms:	HRF; p02; p23; TCTP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)