

Product datasheet for **HP206855**

DNAJC7 Human qPCR Primer Pair (NM_003315)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7266
Forward Sequence:	GGAGTGGACAAGAATGCCTCTG
Reverse Sequence:	AACCTCAGCACTGGCTCCACTA
ACCN:	<u>BC033772</u> , <u>NM_003315</u> , <u>NM_003315.1</u> , <u>NM_003315.2</u> , <u>NM_003315.3</u> , <u>BC033772.1</u> , <u>BC011837</u> , <u>BC011837.2</u> , <u>BC003601</u> , <u>BP378441</u> , <u>BQ445927</u> , <u>BX647209</u> , <u>NM_003315.4</u>
UniProt ID:	<u>Q99615</u>
Synonyms:	DJ11; DJC7; TPR2; TTC2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)